

Test Update UPEP with Reflex to Immunofixation Electrophoresis

Effective Thursday, October 6, 2016, the collection requirements for Urine Protein Electrophoresis with Reflex to Immunofixation Electrophoresis have been updated.

In order to prevent unnecessary delays in testing, and to ensure proper sample collection and communication of collection information, separate EPIC orders have been created for **UPEP for Timed Urine** collection (requires date and times sample collection started and ended, includes 24 hour urines) and **UPEP for Random Urine** collection.

| Test Mnemonic/Description | UPEP w/ Reflex to IFE, Timed Urine | UPEP w/Reflex to IFE, Random Urine |
|---------------------------|--|---|
| EPIC Code | LAB438 | LAB837 |
| Methodology | Urine protein electrophoresis, Immunofixation electrophoresis | Urine protein electrophoresis, Immunofixation electrophoresis |
| Sample Collection | 10.0 mL urine aliquot of a 24 hour urine collection Mandatory requirement of documenting the date and time at the start of sample collection and date and time at the end of sample collection. | 14 mL random urine collection Two, 7 mL urine collection tubes, no additives (clear top) Minimum volume: 7 mL urine |
| Unacceptable Conditions: | Minimum volume: 6 mL urine Specimen not properly identified, incorrect container, insufficient sample volume, no collection date and times noted for start and end of collection. | Specimen not properly identified, incorrect container, insufficient sample volume |
| СРТ | 84156 Protein Timed Urine; 84166 Protein Electrophoresis and if indicated 86335 IFE Urine | 84156 Protein Urine; 84166 Protein Electrophoresis and if indicated 86335 IFE Urine |

Please contact Gregary Bocsi, DO at 303-848-7050 if you have any questions, or visit our website at https://www.uchealth.org/professionals/Pages/Clinical-Laboratory.aspx for additional information.

Issued: November 1, 2016

Clinical Laboratory
Medical Director Edward R. Ashwood, MD